



- 1) Install utility on Server computer at clinic
- 2) Enter main clinic telephone number, click Next




The screenshot shows a window titled "Download AVImark Sync Installer" with a close button (X) in the top right corner. The main content area contains the text "Clinic Phone Number" above a text input field. The input field has a placeholder format "() - ". At the bottom right of the window, there is a blue button labeled "Next >".

- 3) After installation you will receive a download complete window



The screenshot shows the same window titled "Download AVImark Sync Installer" with a close button (X) in the top right corner. The main content area displays the text "Starting setup program..." followed by "Download complete." Below this text is a solid green horizontal progress bar. At the bottom right of the window, there is a button labeled "Next >".

- 4) The AVImark Sync Install window will open. Check the box to accept the license agreement, click Next.



- 5) Browse to the location of the AVImark folder, click Next
- 6) It will install the AVImark Sync Program and you will receive a window confirming a successful installation. Click Finish to complete setup.

